




# Error-corrected sequencing


Chris Miller, PhD  
c.a.miller@wustl.edu




# Sequencing data is noisy




# Sequencing data is noisy




Error bases, zeros removed (n = 17,043)



Error bases, zeros removed (n = 20,334)



Error bases, zeros removed (n = 20,414)



Error bases, zeros removed (n = 16,777)




# Sources of error

- DNA damage during handling/prep (FFPE)
- PCR Amplification steps
- Base calling on the instrument




# Sequencing data is noisy




Errors are a fact of life




# MDS response to decitabine



# MDS response to decitabine




# MDS response to decitabine






# MDS response to decitabine



# MDS response to decitabine



# Error-corrected sequencing

1) Digest genomic DNA.



2) Hybridize the HaloPlex probe library in presence of the Indexing Primer Cassette. Hybridization results in gDNA fragment circularization and incorporation of indexes and Illumina sequencing motifs.



3) Capture target DNA-probe hybrids. Biotinylation of probe DNA allows capture using streptavidin-coated magnetic beads.



4) PCR amplify targeted fragments to produce a sequencing-ready, target-enriched sample.



# Barcoded sequencing


d. Generate read families




e. Create error-corrected consensus-sequence (ECCS)




f. Align ECCSs




# Ultradeep barcoded sequencing



# MDS response to decitabine



# Ultradeep barcoded sequencing



# Key tools for working with error-corrected reads

- analysis-workflows repo
  - pipelines > alignment\_umi\_molecular.cwl
  - pipelines > alignment\_umi\_duplex.cwl
- fgbio tools
  - <https://fulcrumgenomics.github.io/fgbio/>






# Alignment

- fgbio: ExtractUmisFromBam
  - Pulls UMIs out of sequence and into a tag (or tags)
- Align with BWA
- Picard - MergeBamAlignment
  - Takes the aligned bam, and re-attaches those UMI tags to the reads



# Collapsing/Assembly

- fgbio: GroupReadsByUmi
  - re-sorting the bam to make all reads from the same molecule adjacent
- fgbio: CallMolecularConsensusReads
  - collapse/assemble them into one read per molecule
  - Sets some tags with stats - depth/error rate/etc
- fgbio: FilterConsensusReads
  - Filter based on above tags
- fgbio: ClipBam
  - Remove overlaps




# Results

- A bam with very low error levels
- Can be run through standard variant callers
  - Parameters may need to be tweaked!!




# Results

- Further filtering may be necessary
- Consider site-specific error rates
  - Iterative calling



# Applications


- Circulating tumor DNA (ctDNA)



Early detection!

# Problems

- Clonal Hematopoiesis/Somatic mosaicism




# Screening

- High false positive rate - then what?!



# Source identification

- Tissue specific methylation patterns





# Applications

- Tracking residual disease
  - much closer to "prime time"
- Off-target effects of therapies
  - chemo-induced damage
  - off-target CRISPR/gene editing
- Microbial species
- When you need to find "needle in a haystack"



